

HAVER RT Sample Splitters

for dividing sieve material into representative test sieve samples

RT 75 – RT 50 – RT 37,5 – RT 25 made from galvanized sheet steel

Included in delivery: 1 rack
1 divider
3 receivers with grab handles (8 l)

Sample Splitters	Article-No.
RT 75 (6 divisions 75 mm = 3")	540286
RT 50 (8 divisions 50 mm = 2")	540287
RT 37,5 (12 divisions 37,5 mm = 1½")	540288
RT 25 (16 divisions 25 mm = 1")	540290
Component parts	Article-No.
Rack	540854
Receiver	540847
Divider 6 divisions	541624
Divider 8 divisions	540861
Divider 12 divisions	541718
Divider 16 divisions	540878

RT 12,5 – RT 6,3 made from lacquered sheet steel

Included in delivery: 1 rack made from lacquered sheet steel
1 divider made from stainless steel
3 receivers

Sample Splitters with receivers made from tinned sheet (1,75 l)	Article-No.
RT 12,5 (18 divisions 12,5 mm = ½")	540291
RT 6,3 (12 divisions 6,3 mm = ¼")	540292
Sample Splitters with receivers plastic laminated (2 l)	Article-No.
RT 12,5 (18 divisions 12,5 mm = ½")	540293
RT 6,3 (12 divisions 6,3 mm = ¼")	540294
Component parts	Article-No.
Rack	540892
Receiver made from tinned sheet	540885
Receiver plastic laminated	542920
Divider 18 divisions	540900
Divider 12 divisions	541725
Initial Equipment with mounted „dust-protection“	540187
Retrofit „Dust-protection“ kit	540073

HAVER RPT Rotating Sample Reducer

for dividing sieve material into representative test sieve samples

Included in delivery: 1 divider made from aluminium
1 adjustable vibratory feeder made from stainless steel
Sample bottles with screw lid and funnel

Rotating Sample Reducer	Article-No.
RPT 1:8	540991
RPT 1:10	540992
RPT 1:30	540994

TYLER Sample Reducer RX 18

for sample reduction in the ratio of 16:1

Sample Reducer	Article-No.
RX 18	560053

