


- 347** "Assistent" Rotating mixer, CE
with 5 synchronous turning eccentric PVC-rolls
to take out.

The five PVC-rolls are designed as synchronous rotating surface for sample bottles, especially for samples of whole blood for haematological determinations and for the gentle mixing of other cell suspensions.

For large containers the rolls can be taken out, further optional rolls without eccentric turning (only rotating) can be inserted.

The rotating mixer is also suitable for preventing sedimentation of other samples, such as paint, lacquer, oil, cosmetics.

Specification: 230 volts 50 hertz,
 with cable and safety plug
 36 rpm
 length: 430 mm
 width: 180 mm
 height: 95 mm
 length of rolls: 320 mm
 weight: ca. 3,7 kg
 charge: max. 2.0 kg



347

- 348** "Assistent" Rotating mixer, CE
with 5 synchronous turning eccentric PVC-rolls
as no. 347, but with fixed rolls



348

- 348/20** "Assistent" Rotating mixer, CE
as no. 348, but with shorter rolls

Specifications: Length: 332 mm
 Width: 180 mm
 Height: 95 mm
 length of rolls 200 mm
 230 volts 50 hertz, 32 rpm.
 with cable and safety plug